

Server recommendations by number of concurrent workstations.

****Note**** In regards to 'Dual Core' processor technology. ClinicPro is not a SMP (Symmetric Multiprocessing) program. While a 'Dual Core' system has the advantage of allowing multiple programs to run at one time, it will not offer any additional advantage when compared to a computer with a faster single processor. In fact; some 'Dual Core' systems have a slower base speed than the recommended specifications. Using a computer with a slower processor speed will result in reduced software performance.

Categories:

Small Network: - 2 to 5 systems.
Medium Network: - 5 to 10 systems.
Large Network: - 10 to 20 systems.

Note Please contact ClinicPro Support for clarifications on network size, and if you are considering implementing more than 20 computer systems to run the software.

IMPORTANT

ClinicPro Software does not recommend Windows Vista. [Click Here](#) to read why. Nor will ClinicPro run properly on older versions of Windows, such as 98, NT or 2000. Expect limited to no functionality on these versions due to legacy compatibility issues.

ClinicPro Software does not recommend wireless networking. Due to the inherently slower nature of wireless networks (wireless networks run at least 50% slower than wired networks), expect slower access times from ClinicPro on a wireless network. If you chose to implement a wireless network, we recommend you use hardware optimized to give the fastest speeds possible: Please use a minimum standard 802.11g with WPA encryption (802.11n is also recommended for greater security) and a MIMO based router for stronger signal strength.

General Hardware Requirements:

20+ GB Hardrive
SVGA/XGA Video capable of '1024 by 768' pixel resolution
(CD or DVD) -RW for backup purposes
56k Baud Modem for electronic billing purposes

Recommended Specs for a Small Network Server (2 computers):

Processor: - 2.8+ GHz P4 or 2.0+ GHz 'Dual Core'
RAM: - 1.0 gig (400+ MHz FSB)
Network Connectivity: - 100mb Ethernet
Operating System: - Windows XP Pro, Windows 2003 Server

Recommended Specs for a Medium Network Server (3-5 computers):

Processor: - 3.0+ GHz P4/Xeon or 2.6+ GHz 'Dual Core'
RAM: - 1.5 to 2.0 gig (400+ MHz FSB)
Network Connectivity: - 100/1000mb Ethernet
Operating System: - Windows 2003 Server, Windows XP Pro

Recommended Specs for a Large Network Server (5 + machines):

Processor: - 3.0+ GHz 'Dual Core' Xeon
RAM: - 2.0+ gig (400+ MHz FSB)
Network Connectivity: - 100/1000mb Ethernet (recommend gigabit Ethernet)
Operating System: - Windows 2003 Server Standard Ed.

Printer Recommendation:

HP Laser products are our standard for paper billing and document printing. ClinicPro does have the capability to adjust HCFA print alignment and is compatible with most manufacturer windows based printers.

Recommendations for Internet Access:

Clinic Pro now requires internet access for online updating of the program. We highly recommend High Speed Internet Access (Broadband, Cable, DSL, T1, etc) for the best performance. Dial-Up Internet can be used as an alternative, but it does limit the availability of online training and remote support.

Workstation Recommended Hardware:

General Hardware Requirements:

Processor: - 2.8+ GHz P4 or 2.0+ GHz 'Dual Core'
RAM: - 1.0 gig
Network Connectivity: - 100mb Ethernet
Operating System: - Windows XP Pro
20+ GB Hardrive
SVGA/XGA Video capable of '1024 by 768' pixel resolution