


## Customizing Buttons




The Category button was discussed on the previous section. Let's briefly discuss the following three buttons: Heading, Numeric Reading, and Question /Answer.



The heading will appear above the paragraph for this category if values have been added. When you add this control, it will not be visible. To see this control, or edit its text, you have to click on the “Verbiage”



Button. Once you have created the Category, then add the values to the category, you can add the Heading.




This control will allow you to add numeric values on the screen. Two controls are created when you use this button; a label which contains the text for the reading, and the reading itself. When you double-click in the numeric reading control (box), a number pad will be displayed. It would be optimum to first create a Category for Numeric Reading, and then use this numeric reading button to set your label and control. This is not necessary, use can just create a numeric reading control, but using the Category will add an entire category control for the numeric readings.



Two controls are created when you use this button; a label which contains the text for the question, and a control box to type the free text.


The next three buttons we will discuss are: **Option**, **Context**, **Slider**, and **Launch**.



This button allows you to add option buttons under categories. The option buttons are typically “salmon” color and serve to describe or specify a context button. As discussed previously, there are two types of Categories that can be created to contain your options: Inclusive or Exclusive.



This customization feature is an upcoming feature that will be added with the Patient Kiosk. It will allow you to use this slide feature to chart measures.



Context

**Chief Complaint**

Neck pain

Radiating into right arm

Radiating into left

Cervical

Occiput

C1

These are typically the gray buttons on the EMR template. In most cases, the context controls tell you where something is happening. This includes various body areas, or various symptoms.

The “Context” buttons work in tandem with the “Option” buttons. The context controls might list various body areas and the options buttons would report what type of pain is being experienced by the patient. A different scenario would be the context buttons listing various symptoms and the option buttons to modify or provide additional information about those symptoms. Context buttons require the use of option buttons to report on the encounter SOAP note. ClinicPro training technicians will make this functionality clearer to you during your training.

Launch

This launch control allows you to do two things: Start an entire page from scratch, creating a template that is uniquely yours or use the ability to launch from one screen to another for quicker navigation. For example, after recording a daily note, you can launch directly into the adjustment screen and record your adjustment.